

# Technical Datasheet

---

**Article No. 11102**

**pH Buffer refill kit pH 10.01  
(+/- 0.01/20 °C, blue) traceable to SRM from NIST**

**For Laboratory use.**

Parameter	Value
Appearance	10 vials à 30 ml
Density (20 °C)	1.01 – 1.02 g/ml
pH 10 °C	10.15
pH 15 °C	10.07
pH 20 °C	10.01
pH 25 °C	9.94
pH 30 °C	9.88
pH 35 °C	9.83
pH 40 °C	9.78
pH 50 °C	9.71
pH 60 °C	9.62

**General information:**

Open the bottle of pH buffer solution only briefly to remove a small amount for calibration. Never return used pH buffer solution to the storage bottle.

**Storage:**

Protect from heat and direct sunlight. Store in cool, dry conditions in well sealed receptacles.