

Technical Datasheet

Article No. 11103

**pH Buffer refill kit pH 4.01
(+/- 0.01/20 °C, red) traceable to SRM from NIST**

For Laboratory use.

Parameter	Value
Appearance	10 vials à 30 ml
Density (20 °C)	1.01 – 1.02 g/ml
pH 10 °C	4.00
pH 15 °C	4.00
pH 20 °C	4.01
pH 25 °C	4.01
pH 30 °C	4.02
pH 35 °C	4.02
pH 40 °C	4.03
pH 50 °C	4.05
pH 60 °C	4.08

General information:

Open the bottle of pH buffer solution only briefly to remove a small amount for calibration. Never return used pH buffer solution to the storage bottle.

Storage:

Protect from heat and direct sunlight. Store in cool, dry conditions in well sealed receptacles.