

Technical Datasheet

Art.-No.: 1113

Buffer solution pH 3.00 (pH 3.00 ($\pm 0.02/20$ °C))
traceable to SRM from NIST

For laboratory use.

Parameter	Value
Density	1.01 g/ml
Citrate-NaCl-NaOH	
pH (0 °C)	3.03
pH (10 °C)	3.02
pH (20 °C)	3.00
pH (25 °C)	3.00
pH (30 °C)	3.00
pH (40 °C)	2.99
pH (50 °C)	2.98
pH (60 °C)	2.98
pH (70 °C)	2.98
pH (80 °C)	2.98
pH (90 °C)	2.97