

TECHNICAL DATASHEET

Article no. 1119

Buffer solution pH 9.00 (pH 9.00 (±0.02/20 °C)) traceable to SRM from NIST

For laboratory use.

Parameter	Value
H ₃ BO ₃ -KCl-NaOH	
pH (0 °C)	9.24
pH (10 °C)	9.11
pH (20 °C)	9.00
pH (25 °C)	8.97
pH (30 °C)	8.93
pH (40 °C)	8.86
pH (50 °C)	8.80
pH (60 °C)	8.75
pH (70 °C)	8.71
pH (80 °C)	8.67
pH (90 °C)	8.64