

Technical Datasheet

Article No. 1127

Buffer solution pH 7.00 colour coded (green)
pH 7.00 ($\pm 0.02/20$ °C)
standardized against NIST reference material

For laboratory use.

Parameter	Value
Appearance	green liquid
pH	7.00
Melting point	0 °C
Boiling point	100 °C
Density	1.01 g/ml
<hr/>	
$\text{KH}_2\text{PO}_4\text{-Na}_2\text{HPO}_4\text{-NaN}_3$	
pH (19–21 °C)	6.98–7.02