

# Technical Datasheet

---

## Article No. 1161

Technical buffer solution pH 2.00 colour coded (red)  
pH 2.00 ( $\pm 0.02/25$  °C), traceable to SRM from NIST

For laboratory use.

Parameter	Value
Appearance	red liquid
10 °C	pH = 2.01
15 °C	pH = 2.00
20 °C	pH = 2.00
25 °C	pH = 2.00
30 °C	pH = 1.99
35 °C	pH = 1.99
40 °C	pH = 1.98
50 °C	pH = 1.98
60 °C	pH = 1.98