

Technical Datasheet

Article No. 1165

Technical buffer solution pH 4.01 colour coded (red)
pH 4.01 ($\pm 0.02/25$ °C), traceable to SRM from NIST

For laboratory use.

Parameter	Value
Appearance	red liquid
10 °C	pH = 4.00
15 °C	pH = 4.00
20 °C	pH = 4.00
25 °C	pH = 4.01
30 °C	pH = 4.01
35 °C	pH = 4.02
40 °C	pH = 4.03
50 °C	pH = 4.06
60 °C	pH = 4.10