

Technical Datasheet

Article No. 1167

Technical buffer solution pH 7.00 colour coded (yellow)
pH 7.00 ($\pm 0.02/25$ °C), traceable to SRM from NIST

For laboratory use.

Parameter	Value
Appearance	yellow liquid
10 °C	pH = 7.06
15 °C	pH = 7.04
20 °C	pH = 7.02
25 °C	pH = 7.00
30 °C	pH = 6.99
35 °C	pH = 6.98
40 °C	pH = 6.97
50 °C	pH = 6.97
60 °C	pH = 6.98