

Technical Datasheet

Article no. 1182

2-Propanol puriss., DAC, BP, Belg. Pharm. VI, USP
min. 99.8 %

For laboratory use.

Parameter	Value
Appearance	colourless liquid
Melting point	-89 °C
Boiling point	82 °C
Flash point	12 °C
Density	0.78 g/ml
Molar mass	60.10 g/mol
Non-volatile substances	max. 0.002 %
Water	max. 0.2 %
Free acid	max. 0.002 %
Iron (Fe)	max. 0.0001 %
Heavy metals (as Pb)	max. 0.0002 %
Aldehydes and ketones	max. 0.05 %
Benzene (C ₆ H ₆)	max. 0.0002 %
Methanol (CH ₃ OH)	max. 0.1 %