

Technical Datasheet

Art.-No.: 1356

Sodium thiosulphate solution 0.1 mol/l (0.1 N)

$\text{Na}_2\text{S}_2\text{O}_3 \times 5 \text{H}_2\text{O}$ in H_2O

For laboratory use.

Parameter	Value
Appearance	Colourless liquid
Density	1.00 g/ml
Molar mass	248.18 g/mol
Melting point	0 °C
Boiling point	100 °C
Accuracy (20 °C)	±0.2 %
Traceable to NIST	