

Technical Datasheet

Article No. 1483

2-Methyltetrahydrofuran for synthesis, pure
min. 99.9 % C₅H₁₀O

For laboratory use.

Parameter	Value
Appearance	colourless liquid
Melting point	-136 °C
Boiling point	80.2 °C
Density (20 °C)	0.855 g/cm ³
Molar mass	29.02 g/mol
Refractive index (20 °C)	1.404–1.408
Water (KF)	max. 300 mg/kg
Peroxides (as H ₂ O ₂)	max. 100 mg/kg
Stabilized with ionol	150–400 mg/kg