

# Technical Datasheet

---

**Article No. 1602**

**Potassium dichromate solution 0.02 mol/l for COD determination  
5.884 g K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> /l H<sub>2</sub>O (0.12N)**

**For laboratory use.**

| Parameter     | Value         |
|---------------|---------------|
| Appearance    | Orange liquid |
| Density       | 1.00 g/ml     |
| Molar mass    | 294.19 g/mol  |
| Melting Point | 0 °C          |
| Boiling Point | 100 °C        |