

Technical Datasheet

Article No. 1650

**Kjeldahl tablets NAS
(sodium sulphate 1.0 g / selenium 0.05 g)**

For laboratory use.

Parameter	Value
Appearance	Dark grey speckled, round, flat or FBE tablets
Average weight range	1.008 – 1.092 g
Weight uniformity	max. 10 % outside 0.971 – 1.129 g, none outside 0.945 – 1.155 g
Diameter	12 – 14 mm
Friability	max. 5 %
Nitrogen (N)	max. 0.05 %