

Technical Datasheet

Article No. 1651

**Kjeldahl tablets NAC
(sodium sulphate 1.0 g / copper(II) sulphate 0.1 g)**

For laboratory use.

Parameter	Value
Appearance	Blue speckled tablets
Average weight range	1.065 – 1.155 g
Weight uniformity	max. 10 % outside 1.018 – 1.182 g, none outside 0.990 – 1.210 g
Diameter	12 – 14 mm
Friability	max. 5 %
Nitrogen (N)	max. 0.05 %
Copper sulphate x 5 H ₂ O (CuSO ₄ x 5 H ₂ O)	90 – 110 mg/tablet