

# Technical Datasheet

---

**Article No. 1653**

**Kjeldahl tablets KS 5.0 g  
(potassium sulphate / selenium 100 : 1)**

**For laboratory use.**

Parameter	Value
Appearance	Grey, round, flat tablets
Average weight range	4.90 – 5.10 g
Weight uniformity	max. 10 % outside 4.75 – 5.25 g, none outside 4.65 – 5.35 g
Diameter	18 – 20 mm
Friability	max. 5 %
Nitrogen (N)	max. 0.05 %