

Technical Datasheet

Article No. 1656

**Kjeldahl tablets CL
(potassium sulphate 5.0 g/copper(II) sulphate 0.15 g)**

For laboratory use.

Parameter	Value
Appearance	White to pale blue, round, flat tablets, with blue specks
Average weight range	4.893 – 5.407 g
Weight uniformity	max. 10 % outside 4.764 – 5.530 g, none outside 4.635 – 5.665 g
Diameter	18 – 20 mm
Friability	max. 5 %
Nitrogen (N)	max. 0.05 %