

Technical Datasheet

Article No. 1658

**Kjeldahl tablets CM
(potassium sulphate 3.5 g / copper(II) sulphate 0.1 g)**

For laboratory use.

Parameter	Value
Appearance	Blue speckled, round, flat tablets
Average weight range	3.42 – 3.78 g
Weight uniformity	max. 10 % outside 3.333 – 3.870 g, none outside 3.24 – 5.96 g
Diameter	15 – 17 mm
Friability	max. 5 %
Nitrogen (N)	max. 0.05 %
Copper sulphate x 5 H ₂ O (CuSO ₄ x 5 H ₂ O)	90 – 110 mg/tablet