

Technical Datasheet

Article No. 1673

Potassium hydrogen sulfate puriss.
98.0–102.0 % KHSO₄

For laboratory use.

Parameter	Value
Appearance	Colourless to white solid
Density	2.24 g/cm ³
Molar Mass	136.17 g/mol
Melting point	195 °C
Ammonium (NH ₄)	max. 0.002 %
Calcium (Ca)	max. 0.02 %
Iron (Fe)	max. 0.002 %
Heavy metals (as Pb)	max. 0.002 %
Chloride (Cl)	max. 0.002 %