

Technical Datasheet

Article No. 1918

Hexamethyldisilazan z. A. (HMDS)
min. 97,5 % C₆H₁₉NSi₂

For laboratory use.

Parameter	Value
Appearance	Colourless liquid
Density (20 °C)	0.774 g/ml
Molar mass	161.04 g/mol
Melting point	- 82 °C
Boiling point	126 °C
Flash point	15 °C
Identity	complies
Colour (Hazen)	max. 10
Refractive index (20 °C)	1.4060–1.4090