

Technical Datasheet

Article no.: 2250

Acetic anhydride p. a. (min. 98 %)

C₄H₆O₃

For laboratory use.

Parameter	Value
CAS no.	108-24-7
Appearance/condition	Clear liquid
Colour (APHA)	max. 20
Refractive index (20 °C)	1.3881 – 1.3931
Residue on evaporation	max. 30 ppm
Substances reducing KMnO ₄ (ACS)	complies
Chloride (Cl)	max. 5 ppm
Phosphate (PO ₄)	max. 10 ppm
Heavy metals (as Pb)	max. 2 ppm
Sulphate (SO ₄)	max. 5 ppm
Iron (Fe)	max. 5 ppm
Identity	complies

Version no. 2