


Technical Datasheet

Article No. 2269

Ethylene glycol monobutyl ether for synthesis (min. 99.0 %)



For laboratory use.

Parameter	Value
Appearance	Clear, colourless liquid
Density (20 °C)	0.897 – 0.905 g/ml
Molar mass	118.18 g/mol
Boiling point	167 – 172 °C
CAS-No.	111-76-2
Identity	complies
Refractive index (20 °C)	1.4167 – 1.4207
Residue on evaporation	max. 50 ppm
Water (KF)	max. 0.1 %