

# Technical Datasheet

---

**Article No. 2351**

**Diisopropyl ether p. a. (min. 98.5 %)**  
**C<sub>6</sub>H<sub>14</sub>O**

For laboratory use.

Parameter	Value
Appearance	Clear, colourless liquid
Density (20 °C)	0.719 – 0.729 g/ml
Molar mass	102.18 g/mol
Boiling point	66.5 – 69.5 °C
CAS-No.	108-20-3
Identity (IR)	complies
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	max. 5 ppm
Water (KF)	max. 0.1 %
Alcohol miscibility	complies
Chloroform miscibility	complies
Diethyl ether miscibility	complies
Stabilized with BHT	10 ppm