

Technical Datasheet

Article No. 2414

Calcium carbonate p.a., ISO
min. 99.5 % CaCO₃

For laboratory use.

Parameter	Value
Appearance	colourless solid
Density	2.93 g/cm ³
Mol.-Weight	100.09 g/mol
Melting Point	825 °C
Insoluble matter	max. 0.005 %
Barium (Ba)	max. 0.005 %
Copper (Cu)	max. 0.0005 %
Iron (Fe)	max. 0.001 %
Potassium (K)	max. 0.005 %
Magnesium (Mg)	max. 0.05 %
Sodium (Na)	max. 0.2 %
Lead (Pb)	max. 0.0005 %
Strontium (Sr)	max. 0.05 %
Chloride (Cl)	max. 0.005 %
Sulfate (SO ₄)	max. 0.005 %
Nitrogen (N)	max. 0.001 %
Insoluble matter in hydrochloric acid	max. 0.005 %

Test date:

Date

Use before date:

Date

Technical Datasheet

Article No. 2414

Calcium carbonate p.a., ISO
min. 99.5 % CaCO₃

For laboratory use.

Parameter	Value
-----------	-------

This document has been generated automatically and is valid without signature.