

Technical Datasheet

Article No. 2424

Citric acid solution 25 %
24.0-26.0 % C₆H₈O₇

For laboratory use.

Parameter	Value
Appearance	clear, nearly colourless liquid
Density (20 °C)	1.08 – 1.25 g/cm ³
Mol.-Weight	192.13 g/mol
Boiling range	102 – 108 °C
Sulfated ash	max. 0.05 %
Mercury (Hg)	max. 1 mg/kg
Lead (Pb)	max. 0.5 mg/kg
Arsenic (As)	max. 1 mg/kg
Oxalates	max. 100 mg/kg