

Technical Datasheet

Article No. 2455

Cyclohexane for synthesis
min. 99.8 % C₆H₁₂

For laboratory use.

Parameter	Value
Appearance	colourless liquid
Melting point	6 - 7 °C
Boiling point	80 - 82 °C
Density	0.774 - 0.784 g/ml
Mol.-Weight	84.16 g/mol
Identity	complies
Refractive index (20 °C)	1.414 - 1.4314
Water (KF)	max. 150 ppm
Residue on evaporation	max. 50 ppm