

# Technical Datasheet

---

Article No. 2607

N,O-Bis(trimethylsilyl)acetamide p. a. (BSA)  
 $\text{CH}_3\text{C}[\text{NSi}(\text{CH}_3)_3]\text{OSi}(\text{CH}_3)_3$

For laboratory use.

Parameter	Value
Appearance	Clear, yellowish liquid
Boiling point	71–73 °C
Flash point	40 °C
Molar mass	203.43 g/mol
Identity	complies
Density (20 °C)	0.835 g/ml