

Technical Datasheet

Article No. 2613

Aniline p. a. (min. 99.5 % C₆H₇N)

For laboratory use.

Parameter	Value
Appearance	Light brown liquid
Melting point	-6,2 °C
Boiling point	184 °C
Density	1.02 g/ml
Molar mass	93.13 g/mol
Assay	min. 99.5 %
Water	max. 0.1 %
Nitrobenzene (C ₆ H ₅ NO ₂)	max. 0.003 %
Non-volatile substances	max. 0.005 %