

Technical Datasheet

Article No. 2623

Ammonium persulfate cryst.
min. 99.0 % (NH₄)₂S₂O₈

For laboratory use.

Parameter	Value
Appearance	Colourless solid
Density	1.98 g/cm ³
Molar Mass	228.20 g/mol
Freie Acid (as H ₂ SO ₄)	max. 0.1 %
Chloride (Cl)	max. 0.001 %
Heavy metals (as Pb)	max. 0.00005 %
Iron (Fe)	max. 0.001 %
pH (62 g/100 ml)	5–6
Active oxygen	min. 6.90 %
Moisture(105 °C)	max. 0.1 %
free acid	max. 0.1 %