

Technical Datasheet

Article No. 2862

CHEMSOLUTE® Alkaline liquid cleaner, phosphate-free, non-foaming

Parameter	Value
Appearance	Colourless to slightly yellow liquid
Density	1.08 g/ml
Complexing agent	34.5 – 39.5 %
Free Alkali (as NH ₃)	4.3 – 5.3 %
contains Sodium hydroxide	

CHEMSOLUTE® Alkaline Cleaner, phosphate-free is characterized by its outstanding cleaning power with largely universal effect and minimal foam formation. In addition, the excellent water solubility of all components ensures residue-free cleaning of your laboratory utensils.

- Energetic universal cleaning agent for the main wash cycle.
- Cleans even heavily soiled objects.
- Especially suitable for removing starch and protein residues.
- No foam formation, even with strong circulation of the rinsing solution in the machine