

# Technical Datasheet

Article No. 410

Hydrogen peroxide 25 % puriss. (stab.)  
min. 25 % w/w H<sub>2</sub>O<sub>2</sub>

For laboratory use.

Parameter	Value
Appearance	clear colourless liquid
Density	1.09 g/ml
Molar Mass	34.01 g/mol
Non volatile Substances	max. 0.02 %
Free acid	max. 0.02 %
Heavy metals (as Pb)	max. 0.0005 %
Chloride (Cl)	max. 0.005 %
Phosphate (PO <sub>4</sub> )	max. 0.05 %
Arsenic (As)	max. 0.0003 %
Copper (Cu)	max. 0.002 %
Lead (Pb)	max. 0.001 %
Zinc (Zn)	max. 0.002 %