

# TECHNICAL DATASHEET

Article no. 445

Hydrogen peroxide technical (min. 30.0 %)



For laboratory use.

Parameter	Value
CAS no.	7722-84-1
Appearance of solution	clear, colourless
Density (20 °C)	1.110 - 1.115 g/ml
pH	2.0 - 4.0

**Th. Geyer GmbH & Co. KG**



Dornierstr. 4 – 6  
D-71272 Renningen  
Tel.: +49 7159 1637-0  
Fax: +49 7159 1637-710  
renningen@thgeyer.de  
www.thgeyer.de

BW-Bank (Swift/BIC SOLADEST600)  
IBAN DE85600501010002036302  
Postbank Stuttgart (Swift/BIC PBNKDEFFXXX)  
IBAN DE32600100700000020708  
Deutsche Bank (Swift/BIC DEUTDESSXXX)  
IBAN DE06600700700125518100

St.-Nr. 70093/40018 / USt-IdNr. DE147510304  
Amtsgericht Stuttgart / HRA-Nr. 254140  
Persönlich haftende Gesellschafterin:  
Geyer Beteiligungsgesellschaft mbH  
Amtsgericht Stuttgart / HRB-Nr. 252035  
Geschäftsführer: Lutz-Alexander Geyer / Thomas Roth