

TECHNICAL DATASHEET

Article no.: 6396

Sodium sulphate anhydrous puriss. (min. 99.0 %)

Na₂SO₄

For laboratory use.

Parameter	Specification
Assay	99.0 - 100.5 %
Loss on ignition	max. 0.5 %
Free acid	max. 0.05 %
Free alkali	max. 0.04 %
Arsenic (As)	max. 0.0002 %
Calcium (Ca)	max. 0.005 %
Iron (Fe)	max. 0.005 %
Magnesium (Mg)	max. 0.005 %
Zinc (Zn)	max. 0.01 %
Heavy metals (as Pb)	max. 0.001 %
Chloride (Cl)	max. 0.02 %

This document has been generated automatically and is valid without signature.