


## LABSOLUTE® ERLLENMEYER FLASKS

narrow neck, with beaded rim

### Properties:

- Made of borosilicate glass 3.3
- According to ISO 1773
- With beaded rim and labelling field
- Durable, white volume scale
- Dishwasher proof and autoclaveable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C

### Technical drawing / picture:



The technical drawing is just an example. Technical drawings of all sizes are available on request.

Th. Geyer GmbH & Co. KG

**Value table:**

Item no.	V ml	Scale	D ml	Ø mm	h mm
7.690 025	25	15 – 25	10	42	75
7.690 026	50	20 – 30 – 40 – 50	10	51	90
7.690 027	100	40 – 60 – 80 – 100	20	64	105
7.690 028	200	100 – 150 – 200	50	79	135
7.690 029	250	100 – 150 – 200 – 250	50	85	145
7.690 030	300	100 – 150 – 200 – 250 – 300	50	87	160
7.690 031	500	200 – 300 – 400 – 500	100	105	180
7.690 032	1,000	400 – 600 – 800 – 1000	200	131	220
7.690 033	2,000	1000 – 1500 – 2000	500	166	280
7.690 034	3,000	1500 – 2000 – 2500 – 3000	500	187	310
7.690 035	5,000	2000 – 3000 – 4000 – 5000	500	220	365

**Description of the abbreviations in the value table:**

Item no.	Item number
V	Nominal volume of the Erlenmeyer flask in milliliter (ml)
Scale	Scaling that is printed on the Erlenmeyer flask
D	Division of the scale in milliliter (ml)
Ø	Max. diameter of the Erlenmeyer flask at its widest part in millimeter (mm)
h	Max. height of the Erlenmeyer flask in millimeter (mm)

Further information is part of the declaration of conformity.

**Revision 1.0, Stand: 20.07.2017**

**Th. Geyer GmbH & Co. KG**