

SCREW CAPS & POURING RINGS, PBT

The red LABSOLUTE® screw caps and pouring rings made of PBT are the perfect complement to the LABSOLUTE® clear and amber laboratory glass bottles with GL thread. The caps are delivered with a PTFE-coated silicone seal. Evaporation of highly volatile substances is avoided by this additional seal.



Your advantages

- Temperature resistant up to 200 °C
- Sterilizable with hot air at 180 °C
- Caps are delivered with an additional PTFE coated silicone seal

Please do not hesitate to contact your Th. Geyer sales representative or the product management team (pmlabsolute@thgeyer.de), if there are any questions concerning the chemical resistance of the screw caps and/or pouring rings.

Item no.	Description	Thread	Pieces per pack
7 690 393	Screw cap	GL32	1
7 690 068	Screw cap	GL45	1
7 690 394	Pouring ring	GL32	1
7 690 069	Pouring ring	GL45	1