

SCREW CAPS, GL 32 / GL 45, PP

LABSOLUTE® screw caps made of polypropylene (PP) are the perfect complement to the LABSOLUTE® clear and amber laboratory glass bottles with GL 32 or GL 45 thread.

The standard colour of the screw caps is blue, but you can also get them in the fancy colours yellow, green, grey and purple.



Your advantages:

- Autoclaveable at 121 °C
- No additional seal is required because of the integrated lip seal

Please do not hesitate to contact your Th. Geyer sales representative or the product management team (pmlabsolute@thgeyer.de), if there are any questions concerning the chemical resistance of the screw caps and/or pouring rings.

| Item no. | Colour | Thread GL | Pieces per pack |
|-----------|--------|-----------|-----------------|
| 7 690 391 | Blue | 32 | 1 |
| 7 690 066 | Blue | 45 | 1 |
| 7 690 073 | Yellow | 45 | 1 |
| 7 690 074 | Green | 45 | 1 |
| 7 690 075 | Purple | 45 | 1 |
| 7 690 076 | Grey | 45 | 1 |

Th. Geyer GmbH & Co. KG