

LABSOLUTE® TEST TUBES

thin wall thickness

Properties:

- Made of borosilicate glass 4.9 (also known as FIOLAX® made by Schott)
- According to DIN 12395
- Thin wall thickness without graduation
- Dishwasher proof and autoclaveable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C

Picture:



Value table:

Item no.	V ml	S ml	Ø mm	h mm
7.690 080	2	0.5	8	70
7.690 081	6	0.5	10	100
7.690 082	6.5	0.5	12	75
7.690 083	9	0.5	12	100
7.690 084	16	0.5	14	130
7.690 085	20	0.6	16	130
7.690 086	25	0.6	16	160
7.690 087	35	0.6	18	180
7.690 088	45	0.6	20	180

Description of the abbreviations in the value table:

Item no.	Item number
V	Approximate volume of the test tubes
S	Wall thickness of the test tubes
Ø	Max. diameter of the test tubes without beaded rim
h	Max. length of the test tubes

Further information is part of the declaration of conformity.

Revision 1.0, Stand: 22.08.2017

Th. Geyer GmbH & Co. KG