


LABSOLUTE® FILTER GASKETS

grey

Properties:

- Made of natural rubber
- Conical
- Hardness 40° shore A
- Best economic seal between a filter flask and a Buchner funnel

Technical drawing / picture:



Value table:

Item no.	Size	ϕ_a mm	ϕ_b mm	h mm	s mm
7.690 091	1	21	12	18	2.5
7.690 092	2	27	17	20	3
7.690 093	3	33	21	24	3
7.690 094	4	41	27.5	27	4
7.690 095	5	53	33	34	5
7.690 096	6	68	48	35	5.5
7.690 097	7	78	58	35	6
7.690 098	8	89	66	40	6.5

Th. Geyer GmbH & Co. KG

Description of the abbreviations in the value table:

Item no.	Item number
Size	Size of the filter gasket (pressed in the gaskets during production)
\varnothing_a	Upper diameter of the filter gaskets
\varnothing_b	Lower diameter of the filter gaskets
h	Max. height of the filter gaskets
s	Max. wall thickness of the filter gaskets

One set consisting of one piece of each size can be ordered with item no. 7.690 099.

Revision 1.0, Stand: 10.08.2017

Th. Geyer GmbH & Co. KG