


## LABSOLUTE<sup>®</sup> FLAT BOTTOM FLASKS

one neck, with standard ground joint

### Properties:

- Made of borosilicate glass 3.3
- According to DIN EN ISO 4797
- With standard ground joint
- Permanent, white volume scale & labelling field
- Dishwasher proof and autoclavable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C
- Ideal for a consistent heating of liquids

### Technical drawing / picture:



Th. Geyer GmbH & Co. KG

**Value table:**

<b>Item no.</b>	<b>V ml</b>	<b>NS</b>	<b>D mm</b>	<b>H mm</b>	<b>S mm</b>
7.690 110	50	29/32	51	85	1.0
7.690 111	100	29/32	64	103	1.0
7.690 112	250	29/32	85	130	1.2
7.690 113	500	29/32	105	160	1.3
7.690 114	1,000	29/32	131	187	1.5
7.690 115	2,000	29/32	166	230	1.8

**Description of the abbreviations in the value table:**

Item no.	Item number
V	Nominal volume of the flat bottom flask
NS	Nominal size of the joint
D	Max. diameter of the flat bottom flask
H	Max. height of the flat bottom flask
S	Max. thickness of the wall of the flat bottom flask

Further information is part of the declaration of conformity.

**Revision 1.0, Stand: 01.08.2017**

**Th. Geyer GmbH & Co. KG**