


LABSOLUTE[®] ERLLENMEYER FLASKS with standard ground joint

Properties:

- Made of borosilicate glass 3.3
- According to DIN EN ISO 4797
- With standard ground joint
- Permanent, white volume scale & labelling field
- Dishwasher proof and autoclavable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C

Technical drawing / picture:



Th. Geyer GmbH & Co. KG

Value table:

Item no.	V ml	NS	D mm	H mm	S mm
7.690 210	25	14/23	42	70	1.0
7.690 211	50	14/23	51	85	1.0
7.690 213	100	14/23	64	100	1.2
7.690 212	50	29/32	51	85	1.0
7.690 214	100	29/32	64	100	1.2
7.690 215	200	29/32	79	131	1.2
7.690 216	250	29/32	85	140	1.2
7.690 217	300	29/32	87	156	1.3
7.690 218	500	29/32	105	175	1.3
7.690 219	1,000	29/32	131	220	1.5
7.690 220	2,000	29/32	166	270	1.8

Description of the abbreviations in the value table:

Item no.	Item number
V	Nominal volume of the Erlenmeyer flask
NS	Nominal size of the joint
D	Max. diameter of the Erlenmeyer flask
H	Max. height of the Erlenmeyer flask
S	Max. thickness of the wall of the Erlenmeyer flask

Further information is part of the declaration of conformity.

Revision 1.0, Stand: 01.08.2017

Th. Geyer GmbH & Co. KG