

LABSOLUTE[®] BAFFLED FLASKS

Properties:

- Made of borosilicate glass 3.3
- Four bottom baffles
- Straight neck for metal caps or GL45 screw thread
- Erlenmeyer shape
- Dishwasher proof and autoclavable at 121 °C
- Very high resistance against a wide range of chemicals

Pictures:



Notes:

The metal caps must be ordered separately, e.g. from Duran (9 010 482) or Schütt Biotec (9 231 457).

The baffled flasks with GL45 screw thread are supplied incl. normal blue screw cap and pouring ring. Diaphragm screw caps are available as accessories, e.g. from Duran (6 273 369).

Value table:

Item #	V ml	Neck	Ø mm	H mm
7 690 430	100	38 mm	62	117
7 690 431	250	38 mm	82	138
7 690 432	500	38 mm	101	186
7 690 433	1,000	38 mm	129	222
7 690 434	2,000	38 mm	161	282
7 690 441	250	GL45	82	130
7 690 442	500	GL45	101	165
7 690 443	1,000	GL45	129	205
7 690 444	2,000	GL45	161	265

Description of the abbreviations in the value table:

Item #	Item number
V	Nominal volume of the baffled flasks
Neck	Neck type of the baffled flasks
D	Maximum, outer diameter of the baffled flasks at the bottom
H	Maximum height of the baffled flasks

Revision 1.0, Stand: 26.07.2022