

LABSOLUTE® VOLUMETRIC FLASKS

amber glass, with PE stopper

The LABSOLUTE® volumetric flasks with accuracy class A are made from high-quality amber glazed borosilicate glass 3.3 and have a standard ground joint. They are supplied with PE stoppers. The LABSOLUTE® volumetric flasks comply with DIN EN ISO 1042 and are calibrated to 'In' (20 °C). A batch-specific certificate of conformity can be supplied on request by quoting the AN number on the label and the retrace code on the flask.



Item no.	V ml	NS	Tolerance ± ml	Pieces per pack
7 691 070	5	10/19	0,040	1
7 691 071	10	10/19	0,040	1
7 691 072	20	10/19	0,040	1
7 691 073	25	10/19	0,040	1
7 691 074	50	14/23	0,100	1
7 691 075	100	14/23	0,100	1
7 691 076	200	14/23	0,150	1
7 691 077	250	14/23	0,150	1
7 691 078	500	19/26	0,250	1
7 691 079	1,000	24/29	0,400	1
7 691 080	2,000	29/32	0,600	1

Revision 1.0, Stand: 08.09.2017

Th. Geyer GmbH & Co. KG