


## LABSOLUTE® STANDARD GROUND STOPPERS

Hollow glass stoppers, hexagonal, with tip

### Properties:

- Made of borosilicate glass 3.3
- According to DIN 12252 type C
- Nominal size of the ground joint acc. to DIN 12242
- Dishwasher proof and autoclavable at 121 °C
- Very high resistance against a wide range of chemicals
- High temperature resistance up to 500 °C
- Ideal for sealing glass items with standard ground joint

### Technical drawing / picture:



Th. Geyer GmbH & Co. KG

[www.thgeyer.com](http://www.thgeyer.com)

Dornierstr. 4 – 6 · 71272 Renningen · GERMANY · +49 7159 16370

**Value table:**

Item no.	NS
7.691 091	10/19
7.691 092	12/21
7.691 093	14/23
7.691 094	19/26
7.691 095	24/29
7.691 096	29/32

**Description of the abbreviations in the value table:**

Item no.	Item number
NS	Nominal size of the joint

Revision 1.0, Stand: 17.08.2017

Th. Geyer GmbH & Co. KG