

TECHNICAL DATASHEET

Pasteur pipettes 1 ml - 7 691 200

Features:

Description	Pasteur pipettes
Volume	1 ml
Dimension	See technical drawing
Material	Low Density Poly Ethylene (LDPE)
Graduation	0.25 ml, 0.5 ml, 0.75 ml, 1 ml
Drop volume	Ca. 0.04 ml
Sterilization	-
Purpose	Used for liquid transfer
Shelf life	-

Packaging

Sub packaging	500 pcs. in bulk
Pieces per pack	500
Dimensions in mm (W x H x T)	305 x 190 x 160
Weight per Box in kg	Ca. 1.0

Revision 2 - 18.10.2017

TECHNICAL DRAWING

Pasteur pipettes 1 ml - 7 691 200

