

FUNNELS, BOROSILICATE GLASS 3.3

LABSOLUTE® funnels with short stem are made of high quality borosilicate glass 3.3. They can also be used to decant and to filter hot liquids. LABSOLUTE® funnels with diameter from 40 to 100 mm fit to almost every application problem. All funnels are manufactured on the basis of ISO 4798. Also small needs of customers can be satisfied, because LABSOLUTE® funnels are sold piecewise.



Your advantages:

- Temperature resistance up to 500 °C
- Autoclaveable
- Resistant against many chemicals except strong bases, hydrofluoric acid and hot phosphoric acid
- Stem of a 60 mm fits into a standard ground joint NS 10/19

Item no.	Ø mm	Length stem mm	For filter Ø mm	Pieces per pack
7 692 015	40	40	55 – 70	1
7 692 016	50	50	70 – 90	1
7 692 017	60	60	70 – 90	1
7 692 018	70	70	110 – 125	1
7 692 019	80	80	125 – 150	1
7 692 020	100	100	150 – 185	1