

## ANALYTICAL FUNNELS, PP

The LABSOLUTE® analytical funnels made from transparent polypropylene (PP) are available in five different sizes and can be used for almost every filtration and filling of low-viscosity liquids. Standard filter paper fits into the LABSOLUTE® analytical funnels. Ribs on the inside ensure a fast filtration. Ribs on the outside avoid air bubbles between filter paper and funnel's surface.



### Your advantages:

- Autoclaveable 121 °C
- Shatter-proof
- Very economic and safe alternative to analytical glass funnels
- Suitable for the use in contact with food

Please send a message to [pmlabsolute@thgeyer.de](mailto:pmlabsolute@thgeyer.de) or ask your sales representative, if there are any questions about the chemical resistance of the powder funnels.

Item no.	Ø mm	length stem mm	for Ø filter mm
7 692 021	35	52	55
7 692 022	55	61	90
7 692 023	75	80	125
7 692 024	100	99	180
7 692 025	160	151	240

Th. Geyer GmbH & Co. KG