

PIPETTE STAND, PP, 94 POSITIONS

The vertical, rotatable LABSOLUTE® pipette stand made from polypropylene (PP) with 94 positions of different sizes is the best choice for a space-saving, clearly arranged and secure storage of your graduated pipettes up to 25 ml and all kind of volumetric pipettes.



Your advantages:

- ✓ Autoclavable at 121 °C
- ✓ Reinforced base plate for a secure stand in the lab
- ✓ Additional, rotatable plate to protect the pipette tips during storage
- ✓ Fast and easy installation
- ✓ Space-saving storage of your pipettes
- ✓ Dimensions Ø approx. 240 mm, height approx. 465 mm

Item no.	Description	Pieces per pack
7 695 718	Pipette stand, 94 positions	1