

LABSOLUTE® NARROW-MOUTH BOTTLES

Properties:

- Made of LDPE
- Food-safe according to Regulation (EC) No. 10/2011, 1282/2011 and 1183/2012, as well as Regulation (EC) No. 1935/2004
- Temperature resistance as usual for LDPE
- Chemical resistance as usual for LDPE

All narrow-mouth bottles are natural colored and delivered without closure. Different types of closures can be ordered separately.

The bottles are suitable to store or transport liquids. Long-term storage of food is possible at or below room temperature without any problems. Calcification to different temperature is possible under consideration of the time.

Value table:

Item no.	V ml	Ø1 mm	Ø2 mm	GL	H1 mm
7.696 161	20	32	9	14	57
7.696 162	30	35	9	14	66
7.696 163	50	40	13	18	84
7.696 164	100	45	13	18	102
7.696 165	250	60	18	25	138

Description of the abbreviations in the value table:

Item no.	Item no.
V	Nominal volume of the bottle
Ø1	Outer diameter of the bottle
Ø2	Inner diameter of the neck
GL	Thread size
H	Total height of the bottle without closure

Revision 1.0, Stand: 20.03.2018

Th. Geyer GmbH & Co. KG