

POWDER FUNNELS, PP

The unbreakable and very economic LABSOLUTE® powder funnels made from transparent polypropylene (PP) are available in four different sizes and can be autoclaved at 121 °C. The funnels are suitable to transfer powders, granulated materials and highly viscous liquids.



Your advantages:

- Autoclavable at 121 °C
- Parallel stem avoids chocking of the funnel.
- Ribs on the inside lead to a fast filtration.
- Ribs on the outside prevent air blockage during application.

Please send a message to pmlabsolute@thgeyer.de or ask your sales representative, if there are any questions about the chemical resistance of the powder funnels.

Item no.	Ø mm	Ø stem mm	l stem mm	Pieces per pack
7 696 230	65	16	22	1
7 696 231	80	17	29	1
7 696 232	100	22	33	1
7 696 233	150	28	30	1

Th. Geyer GmbH & Co. KG