

LABSOLUTE® WEIGHING BOATS

made from borosilicate glass 3.3

Properties:

- Made from borosilicate glass 3.3 according to ISO 3585
- Funnel shape for a lost-free ejection of the sample
- Dishwasher proof and autoclaveable at 121 °C
- Very high resistance against a wide range of chemicals
- Compatible with all neck sizes from NS 10/19

Picture:



Value table:

Item #	V ml	Ø1 mm	Ø2 mm	L mm
7.696 262	3	8	18	46
7.696 263	6	8	20	60
7.696 264	10	8	22	70

Description of the abbreviations in the value table:

Item #	Item number
V	Approx. volume of the weighing boat
L	Approx. length of the weighing boat without stem
Ø1	Diameter of the stem
Ø2	Approx. body diameter of the weighing boat

Revision 1.0, Stand: 11.05.2022