

TECHNICAL DATASHEET

Centrifugation Tubes 15 ml - 7 696 843

Feature:	
Description	Centrifuge tube with cap, screwed, in racks
Dimension	See technical drawing
Volume	Working volume: 15 ml
Material	Tube: Polypropylene, compliance with USP Class VI Cap: HDPE, compliance with USP Class VI
Color	Tube translucent, black graduation, white writing area Cap: turquoise
Sterilization	Gamma sterilized (SAL 10 ⁻⁶)
Usage	For single use only
Temperature range	Tube: - 80°C to +121°C Cap: - 80°C to +100°C
Autoclavable	Tube: Yes Cap: No
Purity	DNase and RNase Non-pyrogenic
Centrifugation, max RCF	12,000 x g

Packaging

Pieces per rack	25
Pieces per box	500


Revision 1 - 09.02.2022

TECHNICAL DRAWING

Centrifugation Tubes 15 ml

Centrifuge tube 15 ml

Dimensions (Unit:mm)



Valid for 7 696 840, 7 696 841, 7 696 842, 7 696 843 – Revision 1: 09.02.2022